

Critical Care / Detachable Valve- Introducer Kits

Product Code	Description	Units Per Box
--------------	-------------	---------------




PCI Kit/3

496140	4F x 11cm Introducer Without Valve - 21G x 1.0" (2.5cm) Needle - .021"G/W	10
496141	5F x 11cm Introducer Without Valve - 19G x 1.5" (3.8cm) Needle - .028"G/W	10
496142	6F x 11cm Introducer Without Valve - 18G x 2.75" (7cm) Needle - .035" G/W	10
496143	7F x 11cm Introducer Without Valve - 18G x 2.75" (7cm) Needle - .035" G/W	10

Kit Contents:

- | | |
|---|------------------------------------|
| 1- Disposable Towel | 1- 22g x 1.5" (3.8cm) Needle |
| 3- 4" x 4" Gauze Sponges | 1- 5cc Syringe |
| 1- Vessel Dilator | 1- Radiopaque Sheath |
| 1- #11 Mini Scalpel | 1- Guide Wire (FC,SS, 3mmJ and ST) |
| 1- Arterial Needle, Clear Hub (See Listed Size) | |



PCI Kit/2 Standard

496121	5F x 11cm Detachable Hemo Valve And Sheath For .028 G/W	10
496122	6F x 11cm Detachable Hemo Valve And Sheath For .038 G/W	10
496123	7F x 11cm Detachable Hemo Valve And Sheath For .038 G/W	10
496124	8F x 11cm Detachable Hemo Valve And Sheath For .038 G/W	10
496125	9F x 11cm Detachable Hemo Valve And Sheath For .038 G/W	10

PCI Kit/2 Pediatric Length

496126	4F x 6cm Detachable Hemo Valve And Sheath For .021 G/W	10
496127	5F x 6cm Detachable Hemo Valve And Sheath For .028 G/W	10
496128	6F x 6cm Detachable Hemo Valve And Sheath For .038 G/W	10
496129	7F x 6cm Detachable Hemo Valve And Sheath For .038 G/W	10

PCI Kit/2 Feathered Tip Dilator

496308	4F x 11cm Detachable Hemo Valve, Sheath, And Guide Wire*	10
496309	5F x 11cm Detachable Hemo Valve, Sheath, And Guide Wire*	10
496310	6F x 11cm Detachable Hemo Valve, Sheath, And Guide Wire*	10
496311	7F x 11cm Detachable Hemo Valve, Sheath, And Guide Wire*	10
496312	8F x 11cm Detachable Hemo Valve, Sheath, And Guide Wire*	10

* .021" x 80cm (FC, SS, 3mmJ And Straight)