


Specialty Introducers and Accessories - Introducer Kits

Product Code	Description	Units Per Box
--------------	-------------	---------------



191010	45cm Introducer Set With Multipurpose Tip 6F x 45cm One Piece Introducer (Non-Detachable Valve) - For .038" G/W	10
191020	7.5F x 45cm One Piece Introducer (Non-Detachable Valve) - For .038" G/W	10
191030	8F x 45cm One Piece Introducer (Non-Detachable Valve) - For .038" G/W	10



	E-P Introducer, One Piece Style (No Sideport Or Stopcock)	
496455	5F x 11cm One Piece Introducer (Non-Detachable Valve)	10
496456	6F x 11cm One Piece Introducer (Non-Detachable Valve)	10
496457	7F x 11cm One Piece Introducer (Non-Detachable Valve)	10
496458	8F x 11cm One Piece Introducer (Non-Detachable Valve)	10
496476	6F x 11cm Tuohy Borst Introducer (Non-Detachable Valve)	10
496477	7F x 11cm Tuohy Borst Introducer (Non-Detachable Valve)	10
496478	8F x 11cm Tuohy Borst Introducer (Non-Detachable Valve)	10

Kit Contents:

- 1- Introducer
- 1- Guide Wire (.038" x 45cm, FC, SS, 3mmJ & Straight)
- 1- Vessel Dilator



497296	Detachable Hemo Valve (With Stopcock) 4F - 6F Detachable Hemo Valve	10
496297	7F - 9F Detachable Hemo Valve	10



497302	80cm Septishield Catheter Sleeve 80cm Catheter Sleeve To Help Protect Catheter From Touch Contamination When Placed Through An Introducer (For use with All Introducers except CLI Style Introducers With Short Hub, Product codes: 496064-067, 496004-009, 496035-038, 490027-029)	10
--------	--	----

500200	Universal Percutaneous Introducer Tray (Does Not Include Introducer)	10
--------	---	----

Tray Contents:

- | | |
|--------------------------------|---|
| 1- Povidone/Iodine Prep Swab | 1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub |
| 1- Disposable Towel | 2- 5cc Ampule Lidocaine HCL 1% |
| 6 - 4" x 4" Gauze Sponges | 1- .035" x 65cm FC, SS, 3mmJ & ST Guide Wire |
| 1- Fenestrated Drape w/Ripline | 2- 5cc Syringe L/S |
| 1- Injection Site | 1- 8F Vessel Dilator |
| 1- 25g x 5/8" (1.6cm) Needle | 1- 22g x 1.5" (3.8cm) Needle |
| 1- 000 Silk Suture | 1- #11 Mini Scalpel |
| 1- Povidone/Iodine Ointment | 1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter |