





Arterial needles		
Product Code	Description	Units Per Box
		
195312	18G x 2.75" (7 cm) AMC/3 Arterial Needle (For 0.032" - 0.038" Guide Wire)	10
195315	19G x 1.5" (3.8 cm) AMC/3 Arterial Needle (For 0.025 - 0.028" Guide Wire)	10
195311	21G x 1.0" (2.5cm) AMC/3 Arterial Needle (For 0.018" - 0.021" Guide Wire) <i>One-Part Thin Wall Short Bevel Needle</i> <i>With Removable Winged Needle Adapter</i>	10
		
195317	18G x 1.5" (3.8cm) AMC/4 Arterial Needle (For 0.032" - 0.038" Guide Wire)	10
195318	18G x 2.75" (7cm) AMC/4 Arterial Needle (For 0.032" - 0.038" Guide Wire)	10
195319	19G x 1.5" (3.8cm) AMC/4 Arterial Needle (For 0.025" - 0.028" Guide Wire)	10
195355	19G x 2.75 (7cm) AMC/4 Arterial Needle (For 0.035" Guide Wire*)	10
195322	21G x 1.0" (2.5cm) AMC/4 Arterial Needle (For 0.018" - 0.021" Guide Wire)	10
195323	21G x 2.0" (5.1cm) AMC/4 Arterial Needle (For 0.018" - 0.021" Guide Wire)	10
195310	21G x 2.75" (7cm) AMC/4 Arterial Needle (For 0.018" - 0.021" Guide Wire) <i>One-Part Thin Wall Short Bevel Needle</i> <i>With Clear Flash back Hub</i> <i>*Extra Thin Wall Needle</i>	10
		
195352	18G x 2.75" (7cm) Winged Arterial Needle (For 0.032" - 0.038" Guide Wire)	10
195330	19G x 2.75" (7cm) Winged Arterial Needle (For 0.035" Guide Wire*) <i>One-Part Thin Wall Short Bevel Needle</i> <i>With Winged Hub</i> <i>*Extra Thin Wall Needle</i>	10

		Arterial needles	
Product Code	Description		Units Per Box
			
195326	18G x 2.875" (7.3 cm) Percutaneous Entry Needle (For .032" - .038" Guide Wire) <i>Two-Part Thin Wall Tapered Cannula and Inner Arterial Beveled Stylet</i>		10
195351	18G x 2.875" (7.3 cm) Percutaneous Entry Needle With Vented Stylet (For .032" - .038" Guide Wire) <i>Two-Part Thin Wall Tapered Cannula and Vented Arterial Beveled Stylet</i>		10
			
195336	18G x 2.125" (5.4 cm) Cournand Style Needle (For .032" - .038" Guide Wire) <i>Three-Part: Outer Thin Wall Tapered Cannula, Inner Arterial Beveled Needle, and Round-Ended Obuturator</i>		10
			
195350	18Gx2.125" (5.4 cm) Modified Cournand Style Needle (For .032"- .038" GuideWire) <i>Two-Part: Outer Thin Wall Tapered Cannula and Inner Arterial Beveled Needle</i>		10