



BIPOLAR PACING CATHETER

STANDARD PIN



Bipolar Pacing Catheter

Series 600



Product Code	Catheter Size	N0. of Lumens	Usable Length	Introducer Recom'd	Catheter Body Material	Maximum Inflation Capacity	Inflated Balloon Diameter	Electrode Spacing
652/1-90P	5 Fr	2	90 cm	6 Fr	PUR	1.0 cc	8 mm	1.0 cm
652/2-90P	5 Fr	2	90 cm	6 Fr	PUR	1.0 cc	8 mm	2.5 cm
652/1-110P	5 Fr	2	110 cm	6 Fr	PUR	1.0 cc	8 mm	1.0 cm
652/2-110P	5 Fr	2	110 cm	6 Fr	PUR	1.0 cc	8 mm	2.5 cm

Description:

The Bipolar Balloon Pacing Catheter is made of extruded radiopaque polymeric tubing with two lumens and two stainless steel electrodes. Electrode spacing is 1.0 or 2.5 cm apart, depending upon which model you choose.

Packaging:

Each catheter is packaged in a sterile peel-pack tray

Catheters Length Marking:

Each catheter has clear, easy to identify increment markings every 10 cm along catheter body