

Arterial Catheter Kits/Trays - Catheters

Product Code	Description	Units Per Box
--------------	-------------	---------------



20g- Yellow Hub, 18g- Pink Hub

Arterial Catheter Mini-Kit

498101	20g x 6" (15cm) With 21G x 2.0" (5.1cm) Needle And 0.021" (40cm) Guide Wire	10
498105	18g x 6" (15cm) With 19G x 2.25" (5.7cm) Needle And 0.028" (40cm) Guide Wire	10
498107	20g x 3" (7.6cm) With 21G x 1.0" (2.5cm) Needle And 0.021" (25cm) Guide Wire	10
498111	18g x 3" (7.6cm) With 19G x 1.5" (3.8cm) Needle And 0.028" (25cm) Guide Wire	10

Kit Contents:

- 1- Radiopaque PTFE Catheter With Female Luer Hub
- 1- Guide Wire (FC/SS/ST)
- 1- Arterial Needle, Clear Hub

Arterial Catheter Kit

498100	20g x 6" (15cm) With 21G x 2.0" (5.1cm) Needle And 0.021" (40cm) Guide Wire	10
498104	18g x 6" (15cm) With 19G x 2.25" (5.7cm) Needle And 0.028" (40cm) Guide Wire	10
498106	20g x 3" (7.6cm) With 21G x 1.0" (2.5cm) Needle And 0.021" (25cm) Guide Wire	10
498110	18g x 3" (7.6cm) With 19G x 1.5" (3.8cm) Needle And 0.028" (25cm) Guide Wire	10

Kit Contents:

- 1- Disposable Towel
- 1- 5cc Syringe
- 1- #11 Mini Scalpel
- 3- 4" x 4" Gauze Sponges
- 1- 9" (23cm) PVC Extension with Stopcock
- 1- Radiopaque PTFE Catheter With Female Luer Hub
- 1- 000 Silk Suture
- 1- Arterial Needle, Clear Hub
- 1- Guide Wire (FC/SS/ST)

Arterial Catheter Tray


498410B	20g x 6" (15cm) With 21G x 2.0" (5.1cm) Needle And 0.021" (40cm) Guide Wire	10
498412B	18g x 6" (15cm) With 19G x 2.25" (5.7cm) Needle And 0.028" (40cm) Guide Wire	10
498411B	20g x 3" (7.6cm) With 21G x 1.0" (2.5cm) Needle And 0.021" (25cm) Guide Wire	10
498413A	18g x 3" (7.6cm) With 19G x 1.5" (3.8cm) Needle And 0.028" (25cm) Guide Wire	10

Tray Contents:

- 1- Fenestrated Drape 22" x 22" (56cm x 56cm) With Ripline
- 3- 4" x 4" Gauze Sponges
- 4- 2" x 2" Gauze Sponges
- 1- Povidone/Iodine Ointment
- 1- 000 Silk Suture
- 1- 5cc Ampule 1% Lidocaine HCL
- 1- #11 Mini Scalpel
- 1- Radiopaque PTFE Catheter With Female Luer Hub
- 1- Arterial Needle, Clear Hub
- 1- 5cc Syringe
- 1- Guide Wire (FC/SS/3mmJ & ST)
- 1- 25g x 5/8" (1.6cm) Needle

Central Venous Catheter Kits/Trays and TPE Trays - Catheters

Product Code	Description	Units Per Box
--------------	-------------	---------------



16g- Purple Hub, 14g- Olive Hub

Blitt Basic Kit/1 CVC Kit

497681	16g x 8" (20cm) CVC Catheter, Single Lumen	10
498682	14g x 8" (20cm) CVC Catheter, Single Lumen	10

Kit Contents:

1- Disposable Towel	1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter
1- 5cc Syringe	1- Radiopaque Polyurethane CVC Catheter
1- #11 Mini Scalpel	1- 000 Silk Suture
3- 4" x 4" Gauze Sponges	1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub
1- 22g x 1.5" (3.8cm) Needle	1- Guide Wire (0.035" x 45cm FC/SS/3mmJ &Straight)
1- 7F Vessel Dilator	

Blitt Maxi-Kit CVC Kit

497683	16g x 8" (20cm) CVC Catheter, Single Lumen	10
497661	16g x 6" (15cm) CVC Catheter, Single Lumen	10
497688	14g x 8" (20cm) CVC Catheter, Single Lumen	10
497660	14g x 6" (15cm) CVC Catheter, Single Lumen	10

Kit Contents:

1- Fenestrated Drape	2- 5cc Syringe
1- 25g x 5/8" (1.6cm) Needle	1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub
1- Povidone/Iodine Ointment	1- Povidone/Iodine Prep Swab
4- 4" x 4" Gauze Sponges	1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter
1- 22g x 1.5" (3.8cm) Needle	1- Guide Wire (0.035" x 45cm FC/SS/3mmJ &Straight)
1- 7F Vessel Dilator	1- #11 Mini Scalpel
1- 000 Silk Suture	1- Radiopaque Polyurethane CVC Catheter

TPE Catheter Tray

008566	5F x 110cm Balloon Flow - Assisted Bipolar Electrode Catheter Tray	5
008416	5F x 115cm Special Care Bipolar Electrode Catheter Tray	5


Tray Contents:

1- #11 Mini Scalpel	2- Disposable Towels
1- 6F Vessel Dilator	1- Guide Wire (0.035" x 45cm FC/SS/3mmJ &Straight)
1- Fenestrated Drape	1- 18g x 2.75" (7cm) Arterial Needle, Clear Hub
5- 4" x 4" Gauze Sponges	1- Povidone/Iodine Prep Swab
1- 25g x 5/8" (1.6cm) Needle	1- 6F Touhy Borst Introducer
1- 21g x 1.5" (3.8cm) Needle	1- Povidone/Iodine Ointment
2- 10cc Syringe L/L	1- 5cc Ampule Lidocaine HCL 1%
2- 5cc Syringe L/L	1- Septishield Catheter Sleeve
1- Bard TPE Catheter Kit	1- 18g x 2.5" (6.4cm) Radiopaque O-T-N Catheter

Bard is a trademark of C.R. Bard, Inc.

Thermodilution And Pulmonary Artery - Catheters


Product Code	Description	Units Per Box
--------------	-------------	---------------



Thermodilution Catheters, Standard (Heparin Coated Balloon)		
497520	7.5F x 110cm 5-Lumen (0.025" G/W) Safety Wedge/Sterile Sleeve	5
497521	7F x 110cm 4-Lumen (0.025" G/W) Safety Wedge/Sterile Sleeve	5
497526	7F x 110cm 4-Lumen (0.025" G/W) Without Safety Wedge Or Sterile Sleeve	5
497537	6F x 110cm 4-Lumen (0.018" G/W) Safety Wedge/Sterile Sleeve	5
497539	7F x 110cm 4-Lumen (0.025" G/W) Safety Wedge	5

Thermodilution Catheters, S-Curve Tip (Heparin Coated Balloon)		
497534	7F x 110cm 4-Lumen (0.025" G/W) Safety Wedge/Sterile Sleeve	5
497536	7F x 110cm 4-Lumen (0.025" G/W) Without Safety Wedge Or Sterile Sleeve	5

Thermodilution Catheters, Soft Body (Heparin Coated Balloon)		
497531	7.5F x 110cm 5-Lumen (0.025" G/W) Safety Wedge/Sterile Sleeve	5



Pulmonary Artery Catheters, Standard (Heparin Coated Balloon)		
497523	7F x 110cm 2-Lumen (0.025" G/W)	5
497527	6F x 110cm 2-Lumen (0.018" G/W)	5
497540	7F x 110cm 2-Lumen (0.035" G/W)	5