

Optical Fiber Scope Catheter

Extra-fine Flexible Angioscope



6000 pixel resolution-color image
Ultimate fineness
Detection of vulnerable plaque
Examination of blood vessel wall

Optical Fiber Scope Catheter



Optical fiber scope catheter is
Extra-fine Flexible Angioscope

Optical fiber scope catheter boasts
6000 Pixel Resolution

Optical fiber scope catheter is a precision examination instrument of
Blood vessel walls especially the coronary arteries

Optical fiber scope catheter is
Efficient for observing the color of plaque

Sterile (EOG)
Non- pyrogenic
Single use only
Light source and video systems are required.

Product Code	Outside diameter (mm)	Effective length (mm)	Length (mm)
OFS-M-2	0.73	1800	3000

Specification

Field of view: 65°
Direction of view: 0°

Material

Optical Glass

